

Guangzhou Netech Environmental Technology Co., Ltd



Add: Room 108, No.18 Mubei West Street, Tianhe District, Guangzhou China

Tel: (8620) +38246009

Contact person: Wendy

Email: sales04@gznetech.com

Skype: netech008

Service-Tel: +8613352827882

Fax: (8620) +32032106

Cell /Whatsapp: +86 13826080743

Service-Mail: service@gznetech.com

Ail-Web: <http://netech21v.en.alibaba.com/>

Com-Web: <http://www.ozonegeneratorsupplier.com/>

Twelve Towers Oxygen Concentrator



1. Patented products

NT-OT series molecular sieve is an air separation plant which has already obtained domestic patent.

2. Rational construction

Modular construction is the main points of our molecular sieve. Through the different combination of oxygen units, the oxygen output can be from 1 LPM to 160 LPM, even higher.

3. Easy to control

Each unit can be work independently, according to the different oxygen output requirement. And free control for the different unit combinations. So that, it can be suitable used for

the standardization of products. Meanwhile, maintenance will be more convenient.

4. The oxygen output is stable and continuously

5. Low-pressure oxygen production working pressure less than 3 KG.

Specifications

| Item | Unit | NTOX3AL | NTOX5AL | NTOX7AL | NTOX10AL | NTOX15AL |
|----------------------|------|--------------------|-----------|-----------|-----------|-----------|
| Oxygen output | LPM | 3 | 5 | 7.5 | 10 | 15 |
| Oxygen concentration | % | 90%±3% | 90%±3% | 90%±3% | 90%±3% | 90%±3% |
| Air flow rate | LPM | 36-45 | 60-75 | 90-115 | 90-115 | 180-200 |
| Pressure(inlet) | Mpa | 0.14-0.18 | 0.14-0.2 | 0.15-0.21 | 0.18-0.25 | 0.35 |
| Pressure(outlet) | Mpa | 0.06-0.08 | 0.06-0.08 | 0.06-0.08 | 0.06-0.08 | 0.06-0.08 |
| Temperature | °C | indoor temperature | | | | |
| Relative humidity | % | ≤ 70%±10% | | | | |
| Noise | Db | ≤60 | ≤60 | ≤60 | ≤60 | ≤60 |
| Dimension | mm | Φ200×310 | Φ200×430 | Φ200×500 | Φ200×560 | Φ200×880 |
| Net weight | kg | 15.70 | 18.30 | 19.30 | 30.60 | 32.5 |

Noted:

1. Prices are including oxygen generating unit, air compressor, and chiller.
2. NTOX15AL-2 comes with 2 units air compressors, 2 chiller.



Guangzhou Netech Environmental Technology Co., Ltd

Add: Room 108, No.18 Mubei West Street, Tianhe District, Guangzhou China

Tel: (8620) +38246009

Fax: (8620) +32032106

Contact person: Wendy

Cell /Whatsapp: +86 13826080743

Email: sales04@gznetech.com

Service-Mail: service@gznetech.com

Skype: netech008

Ail-Web: <http://netech21v.en.alibaba.com/>

Service-Tel:+8613352827882

Com-Web: <http://www.ozonengeneratorsupplier.com/>

Two Towers Oxygen Concentrator

Descriptions:

1. Simple structure, easy to operate, and stable oxygen output.
2. Material: molecular sieve; each unit can work independently, and oxygen output can reach to 3-10 LPM or higher.
3. With chiller for cooling air, no need additional refrigerant dryer for reducing cost.
4. Oxygen purity can up to 90+5%. Two oxygen cylinders exchange working, auto drain water, nitrogen and other impurities output.
5. Parts with shipping: air compressor, chiller.
6. Oxygen outlet pressure: 0.06-0.08Mpa.



| Item | Unit | NT03L | NT05L | NT08L | NT10L |
|----------------------|------|--------------------|------------|-------------|-------------|
| Oxygen output | LPM | 3 | 5 | 8 | 10 |
| Oxygen concentration | % | 90%±3% | | | |
| Power | W | 280 | 450 | / | 650 |
| Pressure(inlet) | Mpa | 0.18-0.2 | 0.2 | 0.2 | 0.2 |
| Pressure(outlet) | Mpa | 0.06-0.08 | | | |
| Temperature | ℃ | indoor temperature | | | |
| Relative humidity | % | ≤85% | | | |
| Noise | Db | ≤ 55 | ≤60 | ≤70 | ≤70 |
| Dimension | mm | 550×160×80 | 550×180×80 | 610×290×120 | 610×310×120 |
| Weight | kg | 9.25 | 16.70 | 18.50 | 19.10 |

Noted:

Prices are including oxygen generating unit, air compressor, and chiller.

Guangzhou Netech Environmental Technology Co., Ltd



Add: Room 108, No.18 Mubei West Street, Tianhe District, Guangzhou China

Tel: (8620) +38246009

Contact person: Wendy

Email: sales04@gznetech.com

Skype: netech008

Service-Tel: +8613352827882

Fax: (8620) +32032106

Cell /Whatsapp: +86 13826080743

Service-Mail: service@gznetech.com

Ail-Web: <http://netech21v.en.alibaba.com/>

Com-Web: <http://www.ozonegeneratorsupplier.com/>

10-40L/Min Industrial Oxygen Generator



Configuration of Oxygen concentrator

1. Built-in oil free air compressor
2. Built-in refrigerated air dryers
3. High efficiency PSA oxygen generate unit
4. Oxygen concentration: 90%±3%
5. Pressure of oxygen outlet: ≤0.1 MPA
6. Control: Door lock, Amp-meter, Working indicator, Power indicator, On/Off, Oxygen flow meter, Voltmeter, Oxygen outlet, Fuse, Power input

Advantage of Oxygen concentrator

1. Low working temperature of air compressor ,thus extending service life of air compressor
2. Built in Small refrigerated air dryer ,thus oxygen concentrator is compact .
3. Oxygen concentrator feed dry, cool ,clean air ,thus extending service life
4. stable oxygen output ,high oxygen purity

Specifications:

| Model | Unit | NT-010LC | NT-015LC | NT-020LC | NT-030LC | NT-040LC |
|------------------------|------|----------|----------|----------|------------|----------|
| Oxygen flow rate | LPM | 10 | 15 | 20 | 30 | 40 |
| Oxygen concentration | % | 90±3 | 90±3 | 90±3 | 90±3 | 90±3 |
| Dimension | cm | 50×84×95 | | | 67×102×950 | |
| Output oxygen pressure | Mpa | ≤0.1 Mpa | | | | |
| Voltage | V | 220V | | | | |
| Power | W | 1060 | 1800 | 2120 | 3180 | 4240 |

Remark:

1. The power of above machine is the rated power.
2. Model NT-010LC and model NT-020LC has stock, the others need to customize.



Guangzhou Netech Environmental Technology Co., Ltd

Add: Room 108, No.18 Mubei West Street, Tianhe District, Guangzhou China

Tel: (8620) +38246009

Contact person: Wendy

Email: sales04@gznetech.com

Skype: netech008

Service-Tel:+8613352827882

Fax: (8620) +32032106

Cell /Whatsapp: +86 13826080743

Service-Mail: service@gznetech.com

Ail-Web: <http://netech21v.en.alibaba.com/>

Com-Web: <http://www.ozonegeneratorsupplier.com/>

NT-OM— Oxygen Generator



Applications:

Hospital, Grass-roots medical, Health institution, Home medical, Nursing home, Resthome, Mine, Public oxygen bar, Beauty, Family health care

Features:

1. Installed two tower oxygen molecular sieve
2. Professional nebulization efficiency, nebulization treatment for respiratory system during intake oxygen
3. Use with pressure control system for oxygen making and supply in one machine, so that avoid to use with oxygen bag and chemical oxygen making equipment, automatic adjusting oxygen separation cycle and oxygen concentration
4. Failure alarm and without electricity alarm
5. Oxygen concentration $93\pm 3\%$;

Specifications:

| Item | Unit | NT-OM3 | NT-OM5 | NT-OM10 |
|------------------------|------|-------------|--------|-----------|
| Oxygen output | LPM | 3 | 5 | 10 |
| Oxygen concentration | % | $93\pm 3\%$ | | |
| Supply power | V Hz | 220 50 | | |
| Power | Kw | 420 | 420 | 550 |
| Oxygen output pressure | Kpa | 0.04-0.07 | | |
| Noise | Db | ≤ 45 | | ≤ 55 |
| NW | Kg | 22 | 23 | 28 |
| Size | mm | 365×375×600 | | |